

Sensor Data Sheet

Ion-Selective Electrode (ISE) Ammonium

DOC 461 7410-E-1.0-DS

Matrix membrane electrode with fluid contact

Article-No. 461 7410



Properties

- Plastic shaft
- Slope 57 ± 2 mV/p NH₄⁺
- Flexible ammonium selective membrane
- Organic ion exchanger in a special solvent, homogeneously distributed in PVC


Applications

- Water management
- Fertilizer detection

Technical Data											
Measuring principle	Ion-Selective Electrode (ISE)										
Measuring range	10 ⁻⁵ – 1 mol/l 0.2 – 18 000 mg/l										
Response time	< 30 s (in 180 mg/l NH ₄ ⁺ solution)										
Stability	± 0.3 mV (20 min); ± 1 mV (24 h)										
Temperature range	0 °C to 40 °C										
Pressure range	max. 1 bar										
pH range	3 pH – 7 pH										
Interfering ions (See also Selectivity Coefficients)	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">Ca²⁺/NH₄⁺ =</td> <td style="width: 50%;">10 000</td> </tr> <tr> <td>Mg²⁺/NH₄⁺ =</td> <td>= 10 000</td> </tr> <tr> <td>Na⁺/NH₄⁺ =</td> <td>1000</td> </tr> <tr> <td>H⁺/NH₄⁺ =</td> <td>100</td> </tr> <tr> <td>K⁺/NH₄⁺ =</td> <td>10</td> </tr> </table> <p style="text-align: center;">ion/NH₄⁺ relation between interfering and measure ions at 10 % error</p>	Ca ²⁺ /NH ₄ ⁺ =	10 000	Mg ²⁺ /NH ₄ ⁺ =	= 10 000	Na ⁺ /NH ₄ ⁺ =	1000	H ⁺ /NH ₄ ⁺ =	100	K ⁺ /NH ₄ ⁺ =	10
Ca ²⁺ /NH ₄ ⁺ =	10 000										
Mg ²⁺ /NH ₄ ⁺ =	= 10 000										
Na ⁺ /NH ₄ ⁺ =	1000										
H ⁺ /NH ₄ ⁺ =	100										
K ⁺ /NH ₄ ⁺ =	10										
Slope	57 ± 2 mV/p NH ₄ ⁺										
Internal resistance	approx. 1 MΩ										
Cable connection	S7 lab connector										
Shaft material	plastic										
Dimensions	Length: 145mm Diameter: 12mm Immersion depth: 120mm										

Sensor Data Sheet

Typical Slope



Selectivity Coefficients

Nicolsky equation: $E = \text{const} + S \cdot \log(a_i + \sum_j K_{i,j} a_j^{n_i/n_j})$

$$a_i = \gamma_i \cdot c_i$$

Ion (i)	Interfering Ion (j)	log $K_{i,j}$
NH ₄ ⁺	Ca ²⁺	-5.0
	Mg ²⁺	-5.0
	Na ⁺	-4.0
	H ⁺	-3.0
	K ⁺	-2.0

Calibration and Measurement

It is recommended to calibrate the electrode with defined solutions. Whether a single-point calibration or a multi-point calibration is required, depends on the type of measurement.

Storage of the Electrode

Store the electrode in a dry environment.